MAGLUMI 600

Chemiluminescence Immunoassay (CLIA) System

Technical Specifications Non-enzyme involved flash chemiluminescence, long stability of reagents **Principle of Luminescence** ABEI label, stable in acid and alkaline buffers Throughput: up to 180 tests/hour Main Features 24 hours ready to use Time to first result:12 minutes Random, Batch and STAT function **Modes of Operation** Up to 16 sample tubes Continuous loading, STAT available Sample Loading Analyzer auto numbered Refrigerated sample area with independent power supply 4 reagents on board Continuous loading Reagent Loading RFID reading all info of reagents Refrigerated reagent area • Integrated kit, ready-to-use, no pretreatment required Integrated kit including calibrator Using superparamagnetic microbeads **Reagent Features** RFID tag storing all info of reagent RFID tag with built-in master curve2-point calibration to adjust master curve Calibration stability: Max 4 weeks Clot detection Liquid level detection for reagent Optional ratio for sample dilution Other Features Auto dilution for high concentration sample Constant 37°C incubation ■ 15-30°C Operating Temperature Windows 7 **Operation System** 10 inch notebook Connection to LIS via COM or network card Interconnection Direct operation on LIS system to test 88×56×50cm, 56kg **Dimensions Weight**



Shenzhen New Industries Biomedical Engineering Co., Ltd (SNIBE Co., Ltd)

4/F,Wearnes Tech Bldg, Shenzhen Science & Industry Park, Nanshan District, Shenzhen, 518057 CHINA

Tel: +86 755 26501514 Fax: +86 755 26654850

Email: sales@snibe.com Web: www.snibe.com







DISTRIBUTOR:



MAGLUMI 600

Chemiluminescence Immunoassay (CLIA) System



www.snibe.com Version 1.0

Outstanding Technology Power of Maglumi

Chemiluminescence Immunoassay (CLIA) System

CLIA uses two important technologies, one is labeling technology which determines reaction mode; and the other is separation technology which determines the sensitivity, accuracy, and precision of the reagents.

Labeling Technology

Two types of labeling technologies are commonly used. One contains enzyme labels and the other contains non-enzyme small molecule labels. Enzyme label reagents are not so stable and are easily affected by the change of storage conditions. Maglumi system applies ABEI labels. ABEI is a non-enzyme small molecule with special molecular formula to enhance stability in acid and alkaline buffer. As fast chemiluminescence, ABEI chemical reaction with sodium hydroxide (NaOH) and Hyperoxide (H₂O₂) finishes the process in 3

Separation Technology

Maglumi uses Nano Magnetic Microbeads, as separation technology, has been widely used in the field of CLIA. Compared with traditional separation technologies, it has advantages as following:

- It helps in enlarging the reaction area of antigens and antibodies to shorten the reaction time.
- Better and faster capture of the antigens and antibodies to enhance the sensitivity.
- As a liquid separation platform, it mixes thoroughly with the reagents to significantly reduce inter-or intra-assay discrepancies.
- Absorbs the antigens or antibodies through chemical reaction to enhance the accuracy.



TEST MENU

Tumor Markers

Ferritin AFP CEA Total PSA f-PSA CA 125 CA 15-3 CA 19-9 PAP CA 50 CYFRA 21-1 CA 242 CA 72-4 NSE S-100 TPA

Cardiac

PG I

PG II

SCCA

CK-MB Troponin I Myoglobin NT-proBNP Aldosterone Angiotensin I Angiotensin II D-Dimer

Drug Monitoring

Cyclosporine A Tacrolimus,FK506

Inflammation Monitoring

PCT CRP

*Under development

FSH LH HCG/β-HCG PRL Estradiol free Estriol Progesterone Testosterone Anti-TPO PAPP-A Rev T3

Hepatic Fibrosis

Thyroid TSH

T4

T3

FT4

FT3

TG

TGA

TRAb

TMA

HA P III P N-P CIV Laminin Cholyglycine

Glyco Metabolism

C-Peptide Insulin IGF-I ICA IAA Proinsulin GAD65

Prenatal Screening

free β-HCG PAPP-A

free Testosterone

Infectious

Fertility

HBsAg HBsAb IgG HBeAg HBeAb IgG HBcAb IgG HCV IgG *HIV *Syphilis

Bone Metabolism

Intact PTH Calcitonin Osteocalcin 25-OH Vitamin D

Anemia

Vitamin B12 Ferritin FA

TORCH

Toxo IgG Toxo IgM Rubella IgG Rubella IgM CMV IgG CMV IgM HSV-1/2 IgG HSV-2 IgG HSV-1/2 IgM

EBV

EBV EA IgG EBV EA IgA EBV VCA IgG EBV VCA IgM EBV VCA IgA EBV NA IgG

Immunoglobulin

IgM IgE

Kidney Function

β_a-MG Albumin

Others

GH Cortisol **ACTH** DHEA-S