


AUTOMATED ANALYSERS



• SELECTRA **ProXS**

Ref. 6003-200


- Throughput : up to 75 photometric tests/hour
- STAT capability
- 48 disposable cuvettes rotor
- 30 reagent positions / 25 sample positions



• SELECTRA **ProS**

Ref. 6003-500

- Throughput : up to 133 photometric tests/hour
up to 333 tests/hour with ISE
- STAT capability
- Washing station for 48 cuvettes rotor
- Water consumption 1L/hour
- 30 reagent positions / 25 sample positions
- Positive Sample Identification (PSID option)
- ISE module (option)
- Integrated touchscreen PC



• SELECTRA **ProM**

Ref. 6003-404

- Throughput : up to 180 photometric tests/hour
up to 360 ISE tests/hour with ISE
- STAT capability
- Washing station for 48 disposable cuvette rotor
- Water consumption 0.5 L/hour
- 32 reagent positions / 62 sample positions
- Positive Sample Identification (PSID option)
- ISE module (option)
- Integrated touchscreen PC



• SELECTRA **ProXL**

Ref. 6003-702

- Throughput : up to 360 photometric tests/hour
up to 540 ISE Tests/hour with ISE
- STAT capability
- Washing station for 48 disposable cuvettes rotor
- Water consumption 1 L/hour
- 64 reagent positions / 110 sample positions
- Positive Sample Identification (PSID)
- ISE module (option)
- Integrated touchscreen PC

SEMI-AUTOMATED ANALYSER



• MICROLAB **300**

Ref. 6002-310

- Simplicity of operating
- Superior optical design
- On-board Quality Control

WHO WE ARE

ELITechGroup is a leading provider bringing high-value diagnostic solutions to laboratories that are closer to the patient.

ELITechGroup comprises the following product lines:

- Clinical Chemistry Systems
- Molecular Diagnostics
- Microbiology

ELITechGroup has become a significant player in the worldwide IVD market:

- Headquarter in Paris, with strong European and US presences
- Direct subsidiaries in France, Italy, UK, Benelux, Brazil, North America, Australia, New-Zeland, Serbia and Switzerland
- Serves customers in more than 100 countries



OUR CLINICAL CHEMISTRY SYSTEMS

ELITechGroup Clinical Systems, with manufacturing operations in France and in The Netherlands, is a global leader in BenchTop Chemistry analysers.

Our family of scalable systems is designed around the BenchTop laboratory's most important needs.

Robust, trouble-free performance comes from decades of refined designs and continued innovation, in both instrument and reagent technologies.



In Vitro diagnostic medical devices : read carefully information available on labelling and/or on the instructions for use and/or in the user manual



WORLDWIDE OFFICES

Please contact your sales representative for terms, conditions and product availability in your country.



ELITech Clinical Systems SAS (France)

ELITechGroup B.V. (The Netherlands)

Headquarters

Australia
BELUX
Brazil
France
Italy

+33 1 41 45 07 10

1800 815 098
+32 92 820 531
+55 27 3025 1415
+33 4 83 36 10 82
+39 011 97 61 91

New Zealand
Serbia
The Netherlands
UK
United States

0800 555 611
+381 11 2467119
+31 313 430 581
+44 1442 869320
+1 800 453 2725



ELITechGroup

EMPOWERING IVD

www.elitechgroup.com
info@elitechgroup.com



Clinical Systems

REAGENTS & INSTRUMENTS



SELECTRA DEDICATED LINE

BIOCHEMISTRY REAGENTS

REFERENCE	PRODUCT	METHOD	PACK SIZE	TEST N°/KIT*
ALSL-0250	ALANINE AMINOTRANSFERASE (ALT/GPT)	IFCC	8 x 20 mL/8 x 5 mL	528
ALBU-0250	ALBUMIN	BROMOCRESOL GREEN	12 x 20 mL	576
ALPI-0230	ALKALINE PHOSPHATASE (ALP)	IFCC	8 x 8.8 mL/4 x 4.8 mL	328
PASL-0230	ALKALINE PHOSPHATASE (ALP)	DGKC/SCE - DEA	4 x 20 mL/4 x 5 mL	320
AMSL-0230	AMYLASE	IFCC TRACEABILITY - CNPG3	6 x 20 mL	342
ASSL-0250	ASPARTATE AMINOTRANSFERASE (AST/GOT)	IFCC	8 x 20 mL/8 x 5 mL	528
BIDI-0250	BILIRUBIN DIRECT	MODIFIED MALLOY EVELYN	8 x 20 mL/8 x 5 mL	528
BITO-0250	BILIRUBIN TOTAL	MODIFIED MALLOY EVELYN	8 x 20 mL/8 x 5 mL	528
CALA-0250	CALCIUM	ARSENazo III	12 x 20 mL	684
CHLO-0250	CHLORIDE	MERCURIC-THIOCYANATE	12 x 20 mL	684
CHSL-0250	CHOLESTEROL	CHOD - PAP	12 x 20 mL	684
HDLL-0230	CHOLESTEROL HDL DIRECT	DIRECT - SELECTIVE DETERGENT	4 x 21 mL/4 x 7 mL	288
LDLL-0230	CHOLESTEROL LDL DIRECT	DIRECT - SELECTIVE DETERGENT	4 x 21 mL/4 x 7 mL	288
CMSL-0230	CREATINE KINASE - MB (CK-MB)	IMMUNO-INHIBITION- IFCC	4 x 20 mL/1 x 20 mL	320
CKSL-0230	CREATINE KINASE (CK NAC)	IFCC	4 x 20 mL/4 x 5 mL	320
CRSL-0250	CREATININE	ENZYMATIC - PAP	8 x 21 mL/8 x 7 mL	664
GHSL-0250	GLUCOSE	HEXOKINASE	8 x 20 mL/8 x 5 mL	456
GPSL-0250	GLUCOSE	GOD - PAP	12 x 20 mL	684
GISL-0250	Y-GLUTAMYLTRANSFERASE (Y-GT)	IFCC - GLUPA-C	8 x 20 mL/8 x 5 mL	576
FEFE-0230	IRON	FERENE	4 x 14.6 mL/4 x 3.9 mL	200
LLSL-0230	LACTATE DEHYDROGENASE (LDH)	IFCC	4 x 20 mL/4 x 5 mL	320
LPSL-0230	LIPASE	COLORIMETRIC, DGGM	4 x 5.6 mL/4 x 3.6 mL	120
MAGX-0230	MAGNESIUM	XYLIDYL BLUE	8 x 8.8 mL	200
PHOS-0230	PHOSPHORUS	PHOSPHOMOLYBDATE	6 x 20 mL	342
PROB-0250	PROTEIN TOTAL (SERUM)	BIURET	12 x 20 mL	684
PRTU-0250	PROTEIN TOTAL (URINE)	PYROGALLOL RED	12 x 20 mL	816
TGML-0250	TRIGLYCERIDES	GPO - PAP	12 x 20 mL	684
URSL-0250	UREA	UREASE - GLDH	8 x 20 mL/8 x 5 mL	528
AUSL-0250	URIC ACID - BI	URICASE - PAP	8 x 20 mL/8 x 5 mL	512
AUML-0250	URIC ACID - MONO	URICASE - PAP	12 x 20 mL	816

SPECIFIC PROTEINS REAGENTS

REFERENCE	PRODUCT	METHOD	PACK SIZE	TEST N°/KIT*
ICRP-0400	CRP	IMMUNO-TURBIDIMETRY	2 x 25 mL/1 x 5 mL	160
CRPW-0230	CRP Wide Range	LATEX - IMMUNO-TURBIDIMETRY	2 x 24 mL/2 x 8 mL	230
IFRT-0230	FERRITIN	LATEX - IMMUNO-TURBIDIMETRY	2 x 20 mL/2 x 5 mL	176
IHAP-0400	HAPTOGLOBIN	IMMUNO-TURBIDIMETRY	2 x 25 mL/1 x 5 mL	100
HBAC-0240	HbA1c	LATEX - IMMUNO-TURBIDIMETRY	1 x 24 mL/1 x 7.6 mL 1 x 0.4 mL/ 4 x 25 mL	89
IIGA-0400	IGA	IMMUNO-TURBIDIMETRY	2 x 25 mL/1 x 5 mL	80
IIGG-0400	IGG	IMMUNO-TURBIDIMETRY	2 x 25 mL/1 x 5 mL	160
IIGM-0400	IGM	IMMUNO-TURBIDIMETRY	2 x 25 mL/1 x 5 mL	100
IMAL-0400	MICROALBUMIN (URINE)	IMMUNO-TURBIDIMETRY	2 x 25 mL/1 x 5 mL	100
IORO-0400	OROSOMUCOID	IMMUNO-TURBIDIMETRY	2 x 25 mL/1 x 5 mL	100
IPAL-0400	PREALBUMIN	IMMUNO-TURBIDIMETRY	2 x 25 mL/1 x 5 mL	194
IRFA-0230	RHEUMATOID FACTOR	LATEX - IMMUNO-TURBIDIMETRY	2 x 20 mL/2 x 5 mL	144
ITRF-0400	TRANSFERRIN	IMMUNO-TURBIDIMETRY	5 x 25 mL/1 x 5 mL	400

SELECTRA DEDICATED LINE

ISE CALIBRATORS, CONTROLS & SOLUTIONS

REFERENCE	PRODUCT	PACK SIZE
ISCA-0250	ISE CALIBRATORS (2 Levels)	6 x 20 mL/6 x 20 mL
ISCT-0046	ISE CONTROL I	10 x 5 mL
ISCT-0047	ISE CONTROL II	10 x 5 mL
ISRS-0800	ISE REFERENCE SOLUTION	1 x 500 mL
ISDI-0250	ISE DILUENT	12 x 25 mL
ISCC-0280	ISE CLEANER/CONDITIONER	6 x 8 mL/3 x 25 mL

SYSTEM SOLUTIONS

REFERENCE	PRODUCT	PACK SIZE
SLHC-5900-EE	ACID SOLUTION	1 L
SLNA-5900-EE	SYSTEM CLEANING SOLUTION	1 L
SLSY-5900-EE	SYSTEM SOLUTION	1 L

EMPTY VIALS

REFERENCE	PRODUCT	PACK SIZE
CVTL-1010	ECS VIALS/CAPS 10 ML	10 Units
CVTL-2530	ECS VIALS/CAPS 25 ML	30 Units
CVTL-5010	ECS VIALS/CAPS 50 ML	10 Units



GENERIC LINE

BIOCHEMISTRY REAGENTS

REFERENCE	PACK SIZE	METHOD
ALBUMIN		
ALBU-0600	2 x 125 mL	BROMOCRESOL GREEN <i>Standard included</i>
ALBU-0700	4 x 250 mL	BROMOCRESOL GREEN <i>Standard included</i>
ALKALINE PHOSPHATASE (ALP)		
PASL-0420	4 x 50 mL/2 x 26 mL	DGKC/SCE - DEA
ALANINE AMINOTRANSFERASE (ALT/GPT)		
ALSL-0430	4 x 50 mL/2 x 26 mL	IFCC
ALSL-0510	5 x 100 mL/1 x 127 mL	IFCC
AMYLASE		
AMSL-0390	1 x 50 mL	IFCC TRACEABILITY CNPG3
AMSL-0400	6 x 50 mL	IFCC TRACEABILITY CNPG3
ASPARTATE AMINOTRANSFERASE (AST/GOT)		
ASSL-0430	4 x 50 mL/2 x 26 mL	IFCC
ASSL-0510	5 x 100 mL/1 x 127 mL	IFCC
BILIRUBIN DIRECT		
BIDI-0600	2 x 100 mL/1 x 50 mL	MODIFIED MALLOY EVELYN
BILIRUBIN TOTAL		
BITO-0600	2 x 100 mL/1 x 50 mL	MODIFIED MALLOY EVELYN
BILIRUBIN DIRECT & TOTAL		
BITD-0600	1 x 100 mL/1 x 100mL/1 x 50 mL	MODIFIED MALLOY EVELYN
CALCIUM		
CALA-0600	2 x 125 mL	ARSENAZO III <i>Standard included</i>
CHLORIDE		
CHLO-0600	2 x 125 mL	MERCURIC-THIOCYANATE <i>Standard included</i>
CHOLESTEROL		
CHSL-0500	6 x 100 mL	CHOD - PAP
CHSL-0700	4 x 250 mL	CHOD - PAP
CHSL-0507	6 x 100 mL	CHOD - PAP <i>Standard included</i>
CHSL-0707	4 x 250 mL	CHOD - PAP <i>Standard included</i>
CHOLESTEROL HDL DIRECT		
HDLL-0380	1 x 60 mL/1 x 20 mL	DIRECT SELECTIVE DETERGENT
HDLL-0390	3 x 60 mL/3 x 20 mL	DIRECT SELECTIVE DETERGENT
CHOLESTEROL LDL DIRECT		
LDLL-0380	1 x 60 mL/1 x 20 mL	DIRECT SELECTIVE DETERGENT
LDLL-0390	3 x 60 mL/3 x 20 mL	DIRECT SELECTIVE DETERGENT
CHOLINESTERASE *		
CHES-0053	16 x 3 mL/4 x 1.5 mL	BUTYRYL-THIOCHOLINE
CREATINE KINASE (CK NAC)		
CKSL-0430	4 x 50 mL/2 x 26 mL	IFCC
CREATINE KINASE - MB (CK-MB)		
CMSL-0410	2 x 50 mL/1 x 26 mL	IMMUNO-INHIBITION - IFCC
CREATININE		
CRCO-0600	2 x 125 mL	JAFFE <i>Standard included</i>
CRCO-0700	4 x 250 mL	JAFFE <i>Standard included</i>
CREATININE		
CRSL-0630	2 x 100 mL/1 x 70 mL	ENZYMATIC - PAP

BIOCHEMISTRY REAGENTS

REFERENCE	PACK SIZE	METHOD
Y-GLUTAMYLTRANSFERASE (Y-GT)		
GISL-0420	4 x 50 mL/2 x 26 mL	IFCC - GLUPA-C
GLUCOSE		
GHSL-0600	5 x 100 mL/1 x 127 mL	HEXOKINASE
GLUCOSE		
GPSL-0500	6 x 100 mL	GOD - PAP
GPSL-0700	4 x 250 mL	GOD - PAP
GPSL-0507	6 x 100 mL	GOD - PAP <i>Standard included</i>
GPSL-0707	4 x 250 mL	GOD - PAP <i>Standard included</i>
IRON		
FEFE-0600	2 x 100 mL/1 x 50 mL	FERENE <i>Standard included</i>
LACTATE *		
LACT-0100	1 x 103 mL/10 x 10 mL	LO -PAP <i>Standard included</i>
LACTATE DEHYDROGENASE (LDH)		
LLSL-0420	4 x 50 mL/2 x 26 mL	IFCC
MAGNESIUM		
MAGX-0600	2 x 100 mL	XYLIDYL BLUE <i>Standard included</i>
PHOSPHORUS		
PHOS-0600	2 x 125 mL	PHOSPHOMOLYBDATE <i>Standard included</i>
PROTEIN TOTAL (SERUM)		
PROB-0600	2 x 125 mL	BIURET <i>Standard included</i>
PROB-0700	4 x 250 mL	BIURET <i>Standard included</i>
PROTEIN TOTAL (URINE)		
PRTU-0600	2 x 125mL	PYROGALLOL RED <i>Standard included</i>
TRIGLYCERIDES		
TGML-0425	6 x 50 mL	GPO - PAP
TGML-0515	6 x 100 mL	GPO - PAP
TGML-0700	4 x 250 mL	GPO - PAP
TGML-0427	6 x 50 mL	GPO - PAP <i>Standard included</i>
TGML-0517	6 x 100 mL	GPO - PAP <i>Standard included</i>
TGML-0707	4 x 250 mL	GPO - PAP <i>Standard included</i>
UREA		
URSL-0420	4 x 50 mL/2 x 26 mL	UREASE - GLDH
URSL-0500	5 x 100 mL/1 x 127 mL	UREASE - GLDH
URSL-0427	4 x 50 mL/2 x 26 mL	UREASE - GLDH <i>Standard included</i>
URSL-0507	5 x 100 mL/1 x 127 mL	UREASE - GLDH <i>Standard included</i>
URIC ACID		
AUML-0420	6 x 50 mL	URICASE - PAP
AUML-0500	6 x 100 mL	URICASE - PAP
AUML-0427	6 x 50 mL	URICASE - PAP <i>Standard included</i>
AUML-0507	6 x 100 mL	URICASE - PAP <i>Standard included</i>

* Lyophilized/Powder

CALIBRATORS - CONTROLS - STANDARDS

CALIBRATORS (Human serum base)		
REFERENCE	PRODUCT	PACK SIZE
CALI-0550	ELICAL 2 Multiparametric	4 x 3 mL
HDLL-0011	CHOLESTEROL HDL	1 x 1 mL
HDLL-0041	CHOLESTEROL HDL	4 x 1 mL
LDLL-0011	CHOLESTEROL LDL	1 x 1 mL
LDLL-0041	CHOLESTEROL LDL	4 x 1 mL

CONTROLS (Human serum base)		
REFERENCE	PRODUCT	PACK SIZE
CONT-0060	ELITROL I Multiparametric	10 x 5 mL
CONT-0160	ELITROL II Multiparametric	10 x 5 mL
CKMB-0900	CK-MB	4 x 3 mL

STANDARDS (Water base)		
REFERENCE	PRODUCT	PACK SIZE
CHOL-0055	CHOLESTEROL Std. 200 mg/dL	4 x 5 mL
GLUP-0055	GLUCOSE Std. 100 mg/dL	4 x 5 mL
PRTU-0022	PROTEIN Std. 100 mg/dL	4 x 2 mL
TRIG-0055	TRIGLYCERIDES Std. 200 mg/dL	4 x 5 mL
URUV-0055	UREA Std. 50 mg/dL	4 x 5 mL
ACUR-0055	URIC ACID Std. 6 mg/dL	4 x 5 mL

SPECIFIC PROTEINS CALIBRATORS & CONTROLS		
REFERENCE	PRODUCT	PACK SIZE
ICRP-0043	CRP CALIBRATOR SET	5 levels x 1 mL
ICRP-0046	CRP CONTROL I	2 x 1 mL
ICRP-0047	CRP CONTROL II	2 x 1 mL
CRPW-0043	CRP WR CALIBRATOR SET	5 x 2 mL
CRPW-0045	CRP WR CONTROL	2 x 2 mL
IFRT-0042	FERRITIN CALIBRATOR	2 x 3 mL
HBAC-0043	HbA1c CALIBRATOR SET	4 x 0.5 mL
HBAC-0049	HbA1c CONTROL LOW/HIGH	2 x 0.5 mL/2 x 0.5 mL
IMAL-0043	MICROALBUMIN (URINE) CALIBRATOR SET	5 levels x 1 mL
IMAL-0046	MICROALBUMIN (URINE) CONTROL I	2 x 1 mL
IMAL-0047	MICROALBUMIN (URINE) CONTROL II	2 x 1 mL
IPRO-0043	PROTEIN CALIBRATOR SET	5 levels x 1 mL
IRFA-0042	RF CALIBRATOR	2 x 2 mL
IRCT-0046	RHEUMATOLOGY CONTROL I	2 x 1 mL
IRCT-0047	RHEUMATOLOGY CONTROL II	2 x 1 mL